

# ANIOSYME FIRST

CE  
0459

Disinfectant detergent  
foam for soiled  
instrumentation



- Collection of instruments and transport while dry: ensures that adequate humidity is maintained
- Aldehyde-free formulation: no fixation of proteins
- Anti-microbial spectrum: intermediate level

## INDICATIONS

Clean and disinfect medical and surgical instruments:

- During transport in non-humid mode
- Prior to machine cleaning

The product maintains an adequate level of humidity to avoid blood and other organic materials drying.

## CHARACTERISTICS

- pH: neutral (approx. 7)
- Ready-to-use, ergonomic foam gun shape
- No production of aerosols, avoiding any risk of emanation from the product
- Neutral pH: compatible with any type of material (including alloys)
- Chloride-free: non-corrosive for materials
- Tri-enzymatic compound: proteases, lipases, amylases, combined with surface agents
- Tried and tested (Degreasing Cleaning Power DCP)

READY-TO-USE  
PRODUCT

