

- **Cleaning disinfectant for medical device**
- **Controlled foaming effect facilitating wiping**



INDICATIONS: Cleaning and disinfection of medical devices.

COMPOSITION: Quaternary ammonium propionate, guanidinium acetate, n-propanol, non ionic surfactant, fragrance, excipients.

MICROBIOLOGICAL PROPERTIES: Bactericidal: EN 1276, NF T 72-190. Bactericidal in dirty conditions: EN 1276, NF T 72-170. Active against *M. tuberculosis* (T.B.). Yeasticidal: EN 1275, EN 13624, EN 14562. Fungicidal: EN 13624, EN 14562. Active against HIV-1, HBV, Rotavirus, Herpesvirus and BVDV (surrogate of HCV).

INSTRUCTIONS FOR USE: Ready to use product. Direct spraying: spray on surface or spread the solution with a clean linen rub, wipe the surface, wait until dry. Do not rinse.

PACKAGINGS:

- 12 bottles, 750 ml each with foaming sprayer



EXCLUSIVE DISTRIBUTOR (Sole Agent) :

PT. PRIMA ALKESINDO NUSANTARA

Gading Bukit Indah
 Jl. Bukit Gading Raya Blok C 12-15
 Kelapa Gading Pemat, Jakarta 14240
 Telp. : (021) 4513874 - 75 Fax. : (021) 4513876
 E-mail : prima@cbn.net.id
 www.primaalkesindo.com