

# ANIOSPRAY 29

- Aldehyde-free formulation, without colouring nor perfuming substance.
- Active against  $H_2N_1$



**INDICATIONS:** Quick disinfection of clean medical devices in operating theatre, highrisk areas, examining rooms, care units, etc.

**COMPOSITION:** Didecyldimethylammonium chloride, polyhexamethylene biguanide, ethanol, glycolic acid.

**MICROBIOLOGICAL PROPERTIES:** Bactericidal: EN 1040, EN 1276, NF T72-171, MRSA (EN 13727). Active against Mycobacterium tuberculosis (T.B.). Yeasticidal (C. albicans). Fungicidal: EN 1275, T 72-301 (A. niger, A. fumigatus). Active against HIV-1, BVDV (surrogate of HCV), Rotavirus, HBV and Influenza virus ( $H_2N_1$  and  $H_5N_1$ ).

**INSTRUCTIONS FOR USE:** Ready to use solution. Has to be used with a hand sprayer or Spraymatic 10 S. Let dry. Do not rinse.

**PACKAGINGS:**

- 12 bottles, 1 L each with 4 sprayers
- 4 cans, 5 L each



**EXCLUSIVE DISTRIBUTOR (Sole Agent) :**

**PT. PRIMA ALKESINDO NUSANTARA**

Gading Bukit Indah  
 Jl. Bukit Gading Raya Blok C 12-15  
 Kelapa Gading Permai, Jakarta 14240  
 Telp. : (021) 4513874 - 75 Fax. : (021) 4513876  
 E-mail : prima@dcbn.net.id  
 www.primaalkesindo.com