



**INDICATIONS:** Airborne disinfection of surfaces of medical devices after cleaning stage (Terminal disinfection).

**COMPOSITION:** Hydrogen peroxide (50 mg/g), excipients.

**MICROBIOLOGICAL PROPERTIES:** Solution is Bactericidal (EN 1040, EN 1276, EN 13727, EN 13697), tuberculocidal (EN 14348), fungicidal (EN 1650, EN 13624, EN 13697) and active against BVDV (surrogate of HCV), against Herpes virus and PRV (surrogate of HBV), Sporicidal (Bacillus sp, Clostridium sp).

**INSTRUCTIONS FOR USE:** Ready to use product. Has to be used without human presence, Aerosolization with AEROSEPT AF. Consumption dose : 12 ml/m<sup>2</sup>. Contact time: 6 hours. Do not rinse.

**PACKAGINGS:**

• 4 cans, 5 L each



**EXCLUSIVE DISTRIBUTOR (Sole Agent) :**

**PT. PRIMA ALKESINDO NUSANTARA**

Gading Bukit Indah  
 Jl. Bukit Gading Raya Blok C 12-15  
 Kelapa Gading Permai, Jakarta 14240  
 Telp. : (021) 4513874 - 75 Fax. : (021) 4513878  
 E-mail : [primed@cbn.net.id](mailto:primed@cbn.net.id)  
[www.primaalkesindo.com](http://www.primaalkesindo.com)