


ASEPTANIOS TERMINAL HPH



INDICATIONS: Airborne disinfection of surfaces of medical devices after cleaning stage (Terminal disinfection).

COMPOSITION: Formaldehyde, N-3-aminopropyl, N-dodecylpropane 1, 3-diamine, ethanol.

MICROBIOLOGICAL PROPERTIES: Bactericidal and fungicidal at 8 ml/m³ in 120 min: NF T 72-281. Virucidal at 8 ml/m³ in 60 min: NF T 72-180. Active against Mycobacterium tuberculosis (TB), Aspergillus and Legionella: NF T 72-281.

INSTRUCTIONS FOR USE: Ready to use product. Dose: 8 ml/m³. Contact time: 2 hours. To be used without human being. Do not rinse.

PACKAGINGS:
• 4 cans, 5 L each

