




Diagnostics

TURKEY

ODAK Blood Grouping Gel Systems




► OBG Forward Group 1 (Forward Blood Grouping)

Identifies ABO and Rh D groups by using gel method. One card involves 2 tests. Includes control well.

Packing : 12 OBG Forward Group 1 Gel Card, 48 OBG Forward Group 1 Gel Card.

Product No: OBGFG1-12, OBGFG1-48




► OBG Forward Group 2 (Forward Blood Grouping)

Identifies ABO and Rh D groups by using gel method. One card involves 2 tests. Does not include control well, instead includes AB.

Packing : 12 OBG Forward Group 2 Gel Card, 48 OBG Forward Group 2 Gel Card.

Product No: OBGFG2-12, OBGFG2-48




► OBG Forward Group 3 (Forward Blood Grouping + CDE)

Identifies ABO, A1, Rh D VI and CDE by using gel method.

Packing : 12 OBG Forward Group 3 Gel Card, 48 OBG Forward Group 3 Gel Card

Product No: OBGFG3-12, OBGFG3-48



► OBG Forward and Reverse Grouping (Forward + Reverse Blood Grouping)

Classification of blood groups should be based on both forward and reverse grouping. The ABO blood group system is determined directly by testing the red blood cells with Anti - A and Anti - B reagents. Confirmation of the test results is provided by testing the serum with known group A1 and group B red blood cells. OBG Forward and Reverse Grouping test allows the detection of ABO, Rh D, Rh D VI and reverse grouping.

Packing : 12 OBG F&R Gel Card, 48 OBG F&R Gel Card

Product No: OBGFG3-12, OBGFG3-48

ODAK