

Blood Bank Refrigerator

FOUNDED 1865 MANUFACTURER OF REFRIGERATION TECHNOLOGY FOR THE HEALTH SECTOR

Philipp Kirsch GmbH

Okenstraße 103 D-77652 OFFENBURG-GERMANY www.kirsch-medical.com









Self-closing doors

Doors left open accidentally are a thing of the past. (Except BL-100, FROSTER-BL-650).



Ventilation slits

Improve the refrigerating machine's ventilation and are easy to clean. Simplifies the integration of the devices in a clean room or into furniture.



EXCLUSIVE DISTRIBUTOR (Sole Agent)

PT. PRIMA ALKESINDO NUSANTARA Kompleks Perkantoran Gading Bukit Indah Blok C No. 12-15, Jl. Bukit Gading Raya, Kelapa Gading Permai, Jakarta 14240 Telp.: (021) 4513875 Fax.: (021) 4513876 E-mail: corporate@prima-alkesindo.com Website: www.prima-alkesindo.com

BL 176



according to DIN 58371

■ PRO-ACTIVE- Control:

Permanent, proactive monitoring of the performance data and alerts in the case of deviations – so that you can take countermeasures in good time before a fault threatens your chilled goods; world's most accurate temperature control in refrigerated areas thanks to two standard PT-1000 sensors.

- External housing made from galvanised sheet steel (rustproof), with white, anti-scratch powder coating. Length of the plug cable: approx. 2.8 m.
- Adjustable feet to compensate for uneven floors.
- Interior made from smooth aluminium with colourless protective coating.
- Interior consists of two drawers of aluminium on telescopic runners with stops. Stacking height: 23 cm.
- Extra-thick 70-mm insulation, made from high-quality, compression-moulded and environmentally friendly material. Energy saving.
- Self-closing door with easy-to-replace plastic magnetic seal frame, lockable.
- **Door stop** on the right-hand side (by default, see illustration), or the left-hand side, can be retrofitted.

- Forced-air cooling, switches off automatically when you open the door, ensures a uniform temperature and minimizes temperature deviations
- Automatic defrosting with time limit and temperature monitoring.
- **Condensate evaporation** in the refrigerating machine room.
- Key switch protected control panel against tampering. Membrane keyboard with digital temperature display and minimum/maximum temperature memory.
- Warning functions with visual and audible alarm signal in the case of power failure and temperature deviations. In the case of power failure, the monitoring unit remains in operation for approx. 30 hours on battery power. Door open alarm after 60 seconds.
- Alarms can be forwarded using the floating contact (e.g., to a mobile phone with optional GSM-Module or to a control centre).
- Data documentation can be read via USB interface with KIRSCH Datanet software.
- Antifreeze against sub-zero temperatures.
- Statically ventilated refrigerating machine, vibration free, hermetically sealed, energy saving, low-noise, easy to service.

Specification

Capacity	170 litres
Temperature setting	+4 °C
Voltage	220 – 240 V, 50/60 Hz
Power consumption	120 watts
Normal consumption	1.22 kWh/24 h
Admissible ambient temperature	from +10 °C to +38 °C
Heat emission (max.)	186 watts
Overall dimensions (including distance from wall)	w x d x h = 70 x 77 x 122 cm
Interior dimensions	w x d x h = $55 \times 52 \times 60$ cm (usable width 3 cm, usable depth 9 cm, usable height 13 cm less)
Exterior dimensions with door open at 90°	w x d = 70 x 136 cm
Drawer inner dimensions	w x d x h = 50 x 34 x 10 cm
Max. load drawer	50 kg
Weight	Net weight 77 kg, gross weight 87 kg

Optional equipment

- **LED illumination** mounted at the top.
- Equipped with castors.
- Drawer insert (set) with 3 adjustable length and 21 cross dividers per drawer.
- Additional length and cross dividers.
- **Drawer** made of stainless steel instead of aluminium.
- Condensate container can be emptied manually if no condensation is required (e.g., in an operating theatre).
- GSM-Module for alarm text message transmission (e.g., to a mobile phone).
- Temperature documentation:
 - PC-KIT-NET
 - 4-pin PT 100 or PT 1000 temperature sensor, class
 1/3 B, including built-in reference body.
 - Disc-type pen recording thermometer, with 7-day measuring cycle, incorporated, with 100 recording discs and 5 felt tip pens.
 - Mechanical pen-recording thermometer with two waxed paper strips.