

# **Pharmaceutical Refrigerator**


**FOUNDED 1865** MANUFACTURER OF REFRIGERATION TECHNOLOGY FOR THE HEALTH SECTOR

#### Philipp Kirsch GmbH

Okenstraße 103 D-77652 OFFENBURG-GERMANY www.kirsch-medical.com









#### **Fast defrosting**

Fast defrosting significantly reduces defrosting periods with minimal temperature variations.



PT. PRIMA ALKESINDO NUSANTARA
Kompleks Perkantoran Gading Bukit Indah Blok C No. 12-15,
Jl. Bukit Gading Raya, Kelapa Gading Permai, Jakarta 14240
Telp.: (021) 4513875 Fax.: (021) 4513876
E-mail: corporate@prima-alkesindo.com
Website: www.prima-alkesindo.com

## **MED 720**

## **PRO-ACTIVE**



according to DIN 13277

- PRO-ACTIVE-Control: Permanent, proactive monitoring of the performance data and alerts in the case of deviations – so that you can take countermeasures in good time before a fault threatens your chilled goods; world's most accurate temperature control in refrigerated areas thanks to two standard PT-1000 sensors.
- External housing made from galvanised sheet steel (rustproof), with white, anti-scratch powder coating. Length of the plug cable: approx. 3 m.
- Equipped with castors.
- **Adjustable feet** to compensate uneven floors.
- Interior made from smooth aluminium with colourless protective coating. Gibs (every 15 mm) for flexible interior arrangement.
- Interior consists of twelve plastic drawers with stops, on roller runners. There are two adjustable length dividers and eighteen adjustable cross dividers per drawer. Storage surface per drawer = 0.23 m<sup>2</sup>.
- Extra-thick 70-mm insulation, made from high-quality, compression-moulded and environmentally friendly material. Energy saving.
- Self-Closing-Door with easy-to-replace plastic magnetic seal frame, lockable.
- Door hinge on the right-hand side (by default, see illustration), or the left-hand side, can be retrofitted.

- Forced-air cooling, switches off automatically when you open the door, ensures a uniform temperature and minimizes temperature deviations.
- Automatic fast defrosting thanks to time-limited and temperature-monitored reversal of the refrigerant circuit.
- Condensate evaporation in stainless steel heated dish beneath the refrigerator floor.
- Key switch protected control panel against tampering. Membrane keyboard with digital temperature display and minimum/maximum temperature memory.
- Data documentation can be read via USB interface with KIRSCH Datanet software
- Warning functions with visual and audible alarm signal in the case of temperature deviations and power failure. In the case of power failure, the monitoring unit remains in operation for approx. 30 hours on battery power.
- Alarms can be forwarded using the potential-free contact (e.g., to a mobile phone with optional GSM-Module or to a control centre).
- Antifreeze against sub-zero temperatures.
- Ventilation-enforced refrigerating machine, hermetically sealed, energy saving, low noise, easy to service.

### Specification

Capacity	700 litres
Temperature setting	approx. +2 °C to +15 °C
Voltage	220 – 240 V, 50 Hz
Power consumption	250 watts
Normal consumption	1,55 kWh/24 h
Admissible ambient temperature	from +10 °C to +38 °C
Heat emission (max.)	533 watts
Exterior dimensions (including door handle and distance from wall)	w x d x h = 77 x 98 x 195.5 cm
Interior dimensions	$w \times d \times h = 62 \times 77 \times 140 \text{ cm}$ (usable depth 10 cm)
Exterior dimensions with door open at 90°	w x d = 77 x 165 cm
Drawer inner dimensions	w x d x h = 56,3 x 59,5 x 5,6 cm
Max. load drawer	24 kg
Weight	Net weight 164 kg, gross weight 190 kg

### Optional equipment

- Glass door with lock.
- **LED-Illumination** mounted on side wall.
- Adjustable feet to compense uneven floors.
- Decorative door frame for attaching the customer's decorative plates.
- Additional length and cross dividers.
- Shelf w x d = 59.7 x 45 cm on rails or overlays instead of drawer (max. load: 40 kg).
- Condensate container can be emptied manually if no condensation is required (e.g., in an operating theatre).
- 60 Hz cooling machinge upon request
- Vario Support, Storage Shelf on roller runners.
- **GSM-Module** for alarm text message transmission (e.g., to a mobile phone).
- **■** Temperature documentation:
  - PC-KIT-NET (network capable version)
  - 4-pin PT 100 or PT 1000 temperature sensor, class
     1/3 B, including built-in reference body (measurement bottle or refrigerating block).
  - Mechanical pen-recording thermometer with two waxed paper strips.