

# **Laboratory Refrigerator**


### **FOUNDED 1865** MANUFACTURER OF REFRIGERATION **TECHNOLOGY FOR THE HEALTH SECTOR**

### **Philipp Kirsch GmbH**

Okenstraße 103 D-77652 OFFENBURG-GERMANY www.kirsch-medical.com









#### Display at eye-level

The display is positioned at eye-level to make it easier to read and simplify operation.



### Optimized air guide plate

The optimized air guide plate reduces temperature deviations inside the system to a minimum.



## **LABO 468**

## PRO-ACTIVE























#### ■ PRO-ACTIVE- Control:

Permanent, proactive monitoring of the performance data and alerts in the case of deviations – so that you can take countermeasures in good time before a fault threatens your chilled goods; world's most accurate temperature control in refrigerated areas thanks to two standard PT-1000 sensors.

- External housing made from galvanised sheet steel (rustproof), with white, anti-scratch powder coating. Length of the plug cable: approx. 2.8 m.
- Adjustable feet to compensate uneven floors.
- Interior made from smooth aluminium with colourless protective coating. Gibs (every 15 mm) for flexible interior arrangement.
- Interior of LABO 340 PRO-ACTIVE consists of one black drawer and four plastic-coated shelves on overlays. Interior of LABO 468 PRO-ACTIVE consists of four plasticcoated shelves on overlays.
- Extra-thick 55 mm energy saving insulation, made from high-quality, compression-moulded and environmentally friendly material.
- Self-closing door with easy-to-replace plastic magnetic seal frame, lockable.
- Door hinge on the right-hand side (by default, see illustration), or the left-hand side, can be retrofitted.

### Specification

	LABO 468 PRO-ACTIVE
Capacity	460 litres
Temperature setting	approx. +0 °C to +15°C
Voltage	220 – 240 V, 50/60 Hz
Power consumption	170 watts
Normal consumption	1.00 kWh/24 h
Admissible ambient temperature	from +10 °C to +38 °C
Heat emission (max.)	280 watts
Exterior dimensions (including door handle and distance from wall)	w x d x h = 74 x 77,5 x 189 cn
Interior dimensions	w x d x h =* 60 x 57 x 128 cm
Exterior dimensions with door open at 90°	w x d = 74 x 142 cm
Shelf size	w x d = 59 x 50 cm
Drawer inner dimensions	-
Max. load drawer/ shelf	–/40 kg
Weight	Net weight 113 kg, gross weight 127 kg

- Forced-air cooling, switches off automatically when you open the door, ensures a uniform temperature and minimizes temperature deviation.
- Automatic defrosting with time limit and temperature monitoring.
- Condensate evaporates in the refrigerating machine
- Key switch protected control panel against tampering. Membrane keyboard with digital temperature display and minimum/maximum temperature memory.
- Warning functions with visual and audible alarm signal in the case of temperature deviations or other malfunctions. Door open alarm after 60 seconds.
- Alarms can be forwarded using a potential-free contact (e.g., to a mobile phone with optional GSM-Module or to a control centre).
- Power failure alarm (visual and audible alarm), the monitoring unit remains in operation for approx. 30 hours on battery power.
- Data documentation can be read via USB interface with KIRSCH Datanet software.
- Antifreeze against sub-zero temperatures.
- Statically ventilated (LABO 340 PRO-ACTIVE) or ventilation-enforced (LABO 468 PRO-ACTIVE) refrigerating machine, hermetically sealed, energy saving, low-noise, easy to service.

### Optional equipment

- Glass door with lock.
- **LED-Illumination** mounted on side wall.
- Decorative door frame for attaching the customer's decorative plates.
- Door fitting.
- **Equipped with castors.**
- Condensate container can be emptied manually if no condensation is required (e.g., in an operating theatre).
- Plastic drawer on roller runners
- Additional length and cross dividers.
- Additional shelf on rails or overlays.
- Basket on rails (only for LABO 468 PRO-ACTIVE) w x d x h = 60 x 45 x 10 cm.
- Aluminium tray with rails (only for LABO 468 PRO-ACTIVE) w x d x h = 60 x 40 x 2.5 cm.
- Vario Support, Storage Shelf on roller runners.
- GSM-Module for alarm text message transmission.
- Temperature documentation:
  - PC-KIT-NET (network capable version)
  - 4-pin PT 100 or PT 1000 temperature sensor, class 1/3 B, including built-in reference body (measurement bottle or refrigerating block).
  - Mechanical pen-recording thermometer with two waxed paper strips.

<sup>\*</sup>usable depth 5 cm less.