

# **Laboratory Refrigerator**


**FOUNDED 1865** MANUFACTURER OF REFRIGERATION **TECHNOLOGY FOR THE HEALTH SECTOR** 

### **Philipp Kirsch GmbH**

Okenstraße 103 D-77652 OFFENBURG-GERMANY www.kirsch-medical.com










### Low-noise compressor

Reduces noise to 40 dB (corresponds to noise level I; sounds perception: whisper).



### High temperature precision

even when installed powerful and ventilated refrigerator by forced-air cooling.



PRIMASINDO

PT. PRIMA ALKESINDO NUSANTARA

Kompleks Perkantoran Gading Bukit Indah Blok C No. 12-15,
Jl. Bukit Gading Raya, Kelapa Gading Permai, Jakarta 14240
Telp.: (021) 4513875 Fax.: (021) 4513876
E-mail: corporate@prima-alkesindo.com

Website: www.prima-alkesindo.com

# LABO 100 **PRO-ACTIVE**





















#### **■** PRO-ACTIVE- Control:

Permanent, proactive monitoring of the performance data and alerts in the case of deviations – so that you can take countermeasures in good time before a fault threatens your chilled goods; world's most accurate temperature control in refrigerated areas thanks to two standard PT-1000 sensors.

- **External housing** made from galvanised sheet steel (rustproof), with white, anti-scratch powder coating. Length of the plug cable: approx. 2.8 m.
- Adjustable feet to compensate uneven floors.
- Interior made from impact-resistant, white plastic with moulded-in shelves.
- **Interior consists** of two plastic-coated shelves.
- Extra-thick 50 mm energy saving insulation, made from high-quality, compression-moulded and environmentally friendly material.
- Lockable door with easy-to-replace plastic magnetic seal frame.
- Door hinge on the right-hand side (by default, see illustration), or the left-hand side, can be retrofitted.
- Forced-air cooling, switches off automatically when you open the door, ensures a uniform temperature and minimizes temperature deviations.

- Automatic defrosting with time and temperature monitoring.
- Condensate evaporates in the refrigerating machine interior.
- Key switch protected control panel against tampering. Membrane keyboard with digital temperature display and minimum/maximum temperature memory.
- Power failure alarm (visual and audible alarm), the monitoring unit remains in operation for approx. 30 hours on battery power.
- Warning functions with visual and audible alarm signal in the case of temperature deviations or other malfunctions. Door open alarm after 60 seconds.
- Alarms can be forwarded using a potential-free contact (e.g., to a mobile phone with optional GSM-Module or to a control centre).
- Data documentation can be read via USB interface with KIRSCH Datanet software..
- Antifreeze against sub-zero temperatures.
- Ventilation-enforced refrigerating machine, hermetically sealed, energy saving, low-noise, easy to service.
- Low-noise compressor reduces noise to 40 dB (corresponds to noise level I; sound perception: whisper).

# Specification

Capacity	95 litres
Temperature setting	approx. +2 °C to +15 °C
Voltage	220 – 240 V, 50/60 Hz
Power consumption	76 watts
Normal consumption	0.64 kWh/24 h
Admissible ambient temperature	from +10 °C to +38 °C
Heat emission (max.)	184 watts
Overall dimensions (including distance from wall)	w x d x h = 54 x 54 x 82 cm
Interior dimensions	w x d x h = $44$ x $43$ x $47$ cm (usable depth 5 cm, usable width 2 cm less)
Exterior dimensions with door open at 90°	w x d = 54 x 106 cm
Shelf size	w x d = 43 x 33 cm
Max. load shelf	25 kg
Weight	Net weight 44 kg, gross weight 53 kg

## Optional equipment

- Glass door with lock.
- Tabletop (exterior dimensions: 85cm in height).
- Decorative door frame for attaching the customer's decorative plates (exterior dimensions: 55cm in width).
- Door fitting.
- **Equipped with castors.**
- Additional plastic drawer w x d x h = 40 x 32 x 5.6 cm.
- Additional length and cross dividers for drawer.
- Additional shelf.
- Condensate container can be emptied manually if no condensation is required (e.g., in an operating theatre).
- GSM-Module for alarm text message transmission (e.g., to a mobile phone).
- Temperature documentation:
  - PC-KIT-NET (network capable version)
  - 4-pin PT 100 or PT 1000 temperature sensor, class 1/3 B, including built-in reference body (measurement bottle or refrigerating block).
  - Mechanical pen-recording thermometer with two waxed paper strips.