

AEROSEPT AF

- Automated aerosol generator
- Process AEROSEPT AF / ASEPTANIOS HP 50 validated: effectiveness proven
- Precision of diffusion thanks to adjustment of real flowrate and required concentration
- Settings of data in French, English and Spanish
- Starting delayed up to 9H, by scheduling

Airborne
disinfection
process

INDICATION

Self-working aerosol generator, mobile and totally automated for disinfection of surfaces from medical devices previously cleaned.

To be used with the ASEPTANIOS HP 50 disinfectant product

RELIABILITY

Automatic suction of the disinfectant product
Ultra-sonic nozzle made of stainless steel,
Compressor watertight totally made of TEFLON,
Peristaltic pump with constant flowrate
(precision: 1%),
Equipment made of single piece in polyethylene.



Laboratoires

ANIOS

La professionnel de la désinfection