


# BlueLine Bag

# delta T

Gesellschaft Für  
Medizintechnik mbH  
GERMANY



## Type :

- BlueLine Bag 5 L
- BlueLine Bag 10 L
- BlueLine Bag 16 L
- BlueLine Bag 30 L

- Insulated bag for individual short term transportation
- Complete with adjustable shoulder strap for easy carrying
- Colour blue, black nylon
- Extra space for transporting documents



The **SilverBag series**-Products that are designed for the safe, efficient and cost effective transportation of blood, blood products, organs, laboratory specimen, medication and pharmaceuticals, such as vaccine, reagent, drugs, etc.

#### **ICE-Akku** (4°C accumulator 1100 ml)



When used in conjunction with the **ICE-Akku** elements, the **SilverBag** keeps the product at the required temperature with an ambient temperature up to +32

The **ICE-Akku** elements are available in a variety of different temperatures :

**-30°C -18°C +4°C  
+22°C and +37°C**

#### **THE SILVER BAG SERIES**

##### **Silver Bag 5 L**

Can store up to 1 pc **ICE-Akku**  
Inner volume : 5 L, weight : 300 g

##### **Silver Bag 10 L**

Can store up to 4 pcs **ICE-Akku**  
Inner volume : 10 L, weight : 500 g

##### **Silver Bag 16 L**

Can store up to 5 pcs **ICE-Akku**  
Inner volume : 16 L, weight : 930 g

Additional PVC-inliner

##### **Silver Bag 30 L**

Can store up to 8 pcs **ICE-Akku**  
Inner volume : 30 L, weight : 1300 g

Additional PVC-inliner

TYPE	max ICE-Akku	duration of 2°C to 10°C
<b>SilverBag 5 L</b>	<b>1 pc</b>	<b>8 hours</b>
<b>SilverBag 10 L</b>	<b>4 pcs</b>	<b>14 hours</b>
<b>SilverBag 16 L</b>	<b>5 pcs</b>	<b>15 hours</b>
<b>SilverBag 30 L</b>	<b>8 pcs</b>	<b>18 hours</b>

The **SilverBag series** is made of robust nylon and insulation material. The bags are designed as insulation for the **ICE-Akku**. The elements fit exactly in the bag and no space is wasted.

**delta T**

Gesellschaft Für  
Medizintechnik mbH  
GERMANY

SOLE AGENT :



**PT. PRIMA ALKESINDO NUSANTARA**

Jl. Bukit Gading Raya Blok C 12-15, Kelapa Gading Permai  
Jakarta Utara 14240  
Telp. : (021) 4513874 - 75 Fax. : (021) 4513876  
E-mail : primed@cbn.net.id