GEL CARD INCUBATOR

CAT. NO.: JKI

Key Benefits:



PID control and industrial sensors provides ultimate results



Only push the start button. Never needs calibration thanks to integrated sensors.

Results

Only in 15 minutes your cards are ready to centrifuge.

Cost Effective

Low startup cost for small to medium labs who require a incubator for their Gel Cards. ODAK Diagnostics also supply gel cards at a competitive price.

Reader and Centrifuge are also available!









INCUBATOR TECHNICAL SPECS

(

TECHNICAL SPECIFICATIONS

DESCRIPTION	SPECIFICATION
Power Supply:	220 V / 50 – 60 Hz
Power requirement:	160 W
Max. Load:	24 Gel Card
Display:	Industrial seven segment displays or LCD Screen
Time Settings:	Could be changed (Precalibrated for ODAK Gel Cards).
Tempature Settings:	Could be changed (Precalibrated for ODAK Gel Cards).
Ambient temperature:	2 °C − 35 °C
Alarm/Buzzer:	Yes
Max. relative humidity:	75 %
Opening the lid:	Manual
Dimensions:	Height: 105 mm / Depth: 400 mm / Width: 270 mm
Weight (with cover box):	11 kg



•