


## OBG AHG GEL



### OBG AHG Gel (Indirect Coombs)

The test contains 8 microtubes Polyspecific anti-human globulin (AHG) gel matrix .The red blood cells will agglutinate in the presence of antibody directed towards the antigen. Agglutinated cells forming a red button on the surface of the gel or agglutinates dispersed in the gel is a positive result and indicates the presence of the corresponding antibody.

Packing : 12 OBG AHG Gel Card, 48 OBG AHG Gel Card.

Product No: OBGAHG-12, OBGAHG-48

